***Influenza Type A Antigen Test Kit***

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## **I. INTRODUCTION**

Influenza type A viruses can infect avian, porcine, equine and other species including humans. Avian Influenza, also known as Fowl Plague, is a viral disease of domestic and wild birds that is characterized by a full range of responses from almost no signs of the disease to very high mortality. Fifteen serologically distinct hemagglutinin and nine neuraminidase subtypes of Influenza type A virus have been isolated from avian species. Subtypes H5 and H7 are associated with significant to catastrophic losses. Disease signs in poultry range from only a slight decrease in egg production to a highly fatal fulminating infection. Signs of infection may include respiratory problems, edema of the head and face or diarrhea.

## **II. INTENDED USE**

The Flu Detect™ Test Strip is an in-vitro, rapid immunochromatographic assay designed to aid in the qualitative detection of Influenza type A virus in tracheal and cloacal specimen from symptomatic birds or flocks. This assay detects all 15 subtypes of Influenza type A virus.

## **III. TEST PRINCIPLE**

This product is developed based on Rapid Immuno-Migration (RIM) technology. The test strip uses two antibodies that are specific to Influenza type A virus. An anti-Influenza A antibody bound to Influenza A antigen present in the sample forms a complex which migrates along a strip and is captured on a sensitized reaction line by the second antibody. The accumulation of the complex causes the formation of a clearly visible pink/purple band. The presence of a control band, located above the reaction line, ensures that the test was performed correctly.

## IV. SAMPLE COLLECTION

Sample collection procedure is very important.

1. Use provided swabs to collect tracheal, oropharyngeal and/or cloacal specimens from avian species (see VII. GENERAL PRECAUTIONS).
2. Samples obtained early in the course of infection will contain the highest detectable amount of virus.
3. Tracheal or oropharyngeal samples should be taken from behind the tongue and into the trachea or oropharyngeal area (not just from the mouth).
4. Cloacal samples should be taken from within the cloacal area avoiding excess solid fecal material or visible blood.
5. Samples should be extracted using provided extraction buffer.
6. Each sample should be extracted with 8 drops (approximately 250 $\mu$ l) of extraction buffer. Excessive extraction buffer volume may result in a decrease in test sensitivity.
7. After collection, samples may be extracted (see IX. PROCEDURE), or swabs may be removed from the handle (by breaking the handle just above the swab) and stored in a capped test tube for testing at a later date (see V. SAMPLE STORAGE).

## V. SAMPLE STORAGE

If samples will not be tested immediately they should be stored at 4°C for up to 24 hours. For prolonged storage, samples should be kept frozen (-20°C). Do not store samples in a self-defrosting freezer. Avoid multiple freeze-thaw cycles.

## VI. KIT STORAGE

Flu Detect™ test kits should be stored between 2 and 30°C (35-86°F). **DO NOT FREEZE**

## VII. GENERAL PRECAUTIONS

1. Do not use kit components after expiration date.
2. Allow samples to come to room temperature (15 to 25°C) (59 to 77°F) before testing.
3. The vial holding the test strips contains a desiccant and should be kept tightly closed when not in use.
4. Use the test strips immediately (within 10 minutes) of removing from desiccant vial.
5. Test strips should only be handled in the upper, labelled region. Avoid handling the center of the test strip.
6. The test strip should be placed in the test tube vertically.
7. Use a separate swab for each sample. Swabs with wooden handles or containing calcium alginate may interfere with the test and cannot be used.
8. Do not centrifuge samples prior to use.
9. Swabs containing visible blood may partly obscure a weak positive band (due to haemoglobin background).

**Note:** *Prior to use, test and control bands appear yellow. The bands are dyed yellow for quality control purposes. The dye does not interfere with the test results and will wash away while the test is developing.*

## VIII. KIT CONTENTS

- A. 1 vial containing 20 test strips
- B. 1 dropper bottle containing extraction buffer (6ml)
- C. 20 swabs
- D. 20 test tubes
- E. 20 caps for test tubes
- F. 1 test tube rack (for 5 tubes)
- G. Instructions for Use

## IX. PROCEDURE

### 1. Sample preparation

- Place 8 drops (approximately 250 $\mu$ L) of extraction buffer into a test tube. (1)
- Swab the trachea or the cloaca of the bird (see IV. SAMPLE COLLECTION).
- Put swab to be tested into collection tube and rotate the swab 5-10 times in the buffer. (2)
- When removing the swab from the tube, press the swab against the side of the tube repeatedly until no more liquid comes from the swab. (3)
- Discard the swab in an appropriate biohazard container.

### 2. Testing Samples

- Remove a test strip from the desiccant vial for each sample to be tested. Handle the test strip on the labelled portion of the strip. (4)
- Place the test strip directly into the test tube containing the sample. Place test strip so that the pink pad is just submerged into the extracted sample. Incubate the test strip in the sample for 15 minutes. (5)
- Remove the test strip from the test tube to read.

**Note:** *If the test strip remains in the test tube for more than 20 minutes a false positive ghost band could appear in place of the reaction band (T).*

### 3. Reading Results

- After 15 minutes, observe the presence or absence of pink/purple bands in the center of the test strip between the two absorption pads.
- The control band appears in the upper end of the test strip (closest to the handle), while the sample test results are read in the lower part of the test strip. (6)
- Discard the test strip in an appropriate biohazard container.

**Note:** *The Control Line on the upper part of the test strip may appear earlier. This does not mean that the test is complete. The test strip must incubate for a full 15 minutes before a sample is interpreted as Negative. The test can be considered to be complete if the Test Line on the lower part of the stick appears before the 15 minute incubation period is over. This sample is interpreted as Positive.*

## **4. Results**

### **Validation**

- The test is validated if the Control Line (pink/purple band) develops in the upper part of the test strip (C). The absence of the Control Line indicates that the test is invalid and must be repeated.

### **Interpretation**

- Negative for Influenza A virus: A single pink/purple band (Control Line) is present in the upper part of the test strip (C).
- Positive for Influenza A virus: Two pink/purple bands (Control Line and Test Line) are present on the test strip (C & T).

### **Positive samples**

- Positive results should be submitted to a reference lab for confirmation and subtype determination.

### **Negative samples**

- Negative results indicate that no detectable Influenza type A virus is present.

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